SWING-BAND® SB3 / SWING-BAND® SB4

FR	2
EN	10
DE	16
EL	23
ES	30
HU	36
IT	42
LV	48
NL	54
PL	61
PT	68
RO	75
AR	82
RU	96
SQ	104

SWING-BAND® SB3 / SWING-BAND® SB4

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
SB IFU
Rev 01 25NOV2022
P 1 / 108

SWING-BAND® SB3 / SWING-BAND® SB4

DEVICES FOR STRESS URINARY INCONTINENCE TREATMENT IN WOMEN

IMPORTANT

- This instruction for use is intended to describe the product proper use. It does not constitute a reference in terms of surgical techniques.
- This device has been designed, tested and manufactured for single patient use. Its reuse, reprocessing and / or re-sterilization can cause malfunctions, injuries, risk of contamination or infection in the patient.

SINGLE USE STERILE PRODUCT
REUSE, REPROCESS OR RESTERILIZE THIS DEVICE IS FORBIDDEN

1 Description

1.1 Intended purpose

SWING-BAND® is a suburethral reinforcement device for stress urinary incontinence treatment. Sterile single-use device consisting of knitted polypropylene sling supplied with or without instrument(s) intended for its placement through the vagina.

1.2 Description

The slings are available in ultra-light SWING-BAND® SB3 (28 g/m²) or light SWING-BAND® SB4 (55 g/m²) versions. Each end of the sling is extended by a traction thread ending in a loop intended to be placed in the eye of the instrument (Figure 1).

A sheath protects the assemblies and ease sling sliding during placement.


Figure 1 : SWING-BAND®

The sling is placed in a tube closed by two caps to avoid soiling the sling during vaginal incisions.

The single-use instrumentation supplied with the implant facilitates the passage of the sling through tissues and allows its placement by the TVT, TVT-O or TOT (Table 1).

The introducer guide guides the passage of the ITG and ITD instruments to the obturator membrane.

Table 1 : Instrument(s) provided as per surgical approach

Range	References	Technique	Pictograms
SWING-BAND® SB3	SB3 10400	Transobturator (TOT, TVT-O) Retropubic (TVT)	
	SB3 10400 SP	Retropubic (TVT)	
	SB3 10400 TO	Transobturator out/in (TOT)	

Range	References	Technique	Pictograms
	SB3 10400 DG	Transobturator in/out (TVT-O)	
SWING-BAND® SB4	SB4 10400	Transobturator (TOT, TVT-O) Retropubic (TVT)	
	SB4 10400 SP	Retropubic (TVT)	
	SB4 10400 TO	Transobturator out/in (TOT)	
	SB4 10400 DG	Transobturator in/out (TVT-O)	

2 Composition

- Sling: Implantable mesh 100% Polypropylene (PP)
Polyethylene Terephthalate (PET) Traction Yarn
Polyolefin sheath
- Tube: Polyethylene Glycol Terephthalic Tube (PET G)
Polyethylene (PE) tube caps
- Instruments: Silicone needle protection (SI)
Instrument handle Polypropylene (PP) / 30% talc
Instrument needles Stainless steel (SS)
Inductor Guide Stainless Steel (SS)

3 Indications and intended population

SWING-BAND® is indicated for adult women with stress urinary incontinence (SUI), pure or mixed, associated with predominant urethral hypermobility.


- After failure of conservative treatments or
- As first-line treatment in cases of severe stress urinary incontinence or
- As second line treatment for recurrent stress urinary incontinence (SUI) after failure of previous surgical treatment.

SWING-BAND® SB3 and SWING-BAND® SB4 differ non significantly by their linear mass. There is no specific indication for one or the other range of slings. The choice of the sling as well as the approach is at the discretion of the surgeon.

4 Contraindications and contraindicated population

SWING-BAND® is contraindicated for growing, pregnant and nursing women or presenting health conditions such as:

- Acute or chronic, local or systemic infections
- Severe muscular, neurological or vascular deficiencies which affect the site concerned
- Any concomitant condition likely to affect the function of the implant
- Allergy or reaction to polypropylene or other component of the product
- Immunodeficiency

5 Performance

5.1 Clinical benefit

SWING-BAND® sling allows total or partial correction of stress urinary incontinence as measured by the patient quality of life improvement. Single-use instruments are used as an aid in combination with SWING-BAND® during surgery, they do not provide direct clinical benefit or performance.

5.2 Clinical performance

SWING-BAND® improves symptoms of the stress urinary incontinence with a complete objective cure rate >92 % and a low stress urinary recurrence rate < 3%.

5.3 Summary of the safety and clinical performance

The summary of the safety and clinical performance characteristics is available on request from the manufacturer via the contact form on its website www.thtbio-science.com. When the European database EUDAMED is functional, it will be available on the following link online, on the EUDAMED database (<https://ec.europa.eu/tools/eudamed>) and accessible by the Basic UDI-DI

Range	Basic UDI-ID
SWING-BAND® SB3	37004863SB3PHIUEF4C
SWING-BAND® SB4	37004863SB4PH2IUEFGK

6 Warnings and precautions for use

The device must be stored in its original sealed packaging away from moisture and light. Check the integrity of the packaging before use. The device should not be used if the packaging has been damaged as its sterility can no longer be guaranteed.

Perform a visual examination of each implant before use, in order to detect any possible defect and discard the product if necessary.

Any incident occurring during the use of the device must be reported to the manufacturer.

Any serious adverse event (SAE) related to the device must be notified to the manufacturer and competent authority of the Member State where the user and / or patient are located.

7 Instructions for use and prerequisites before use

The surgeon must have the appropriate training and experience in urological or gynecological surgery.

SWING-BAND® slings must be used with the appropriate instruments either supplied with the device or available separately.

a. Device opening

The packaging contains: the device, 4 traceability labels and a pre-filled implant card.

Note: Traceability labels has essential information related to the device identification. They are intended for the patient medical file record.

Remove the device, labels and implant cards from the cardboard box. The device is packaged in two peelable pouches which ensure its aseptic presentation up to the sterile field of the operating room. Always open the pouches by the chevron:

- the first one is opened outside the sterile field
- the second one at the sterile field

b. Patient preparation

- The surgery can be performed under local, regional or general anesthesia.
- Place the patient in a dorsal gynecological position, hips flexed, buttocks aligned on the edge of the table.
- Insert a Foley catheter into the urethra to ensure the bladder emptiness.

c. Surgical procedure

Retro-pubic approach (Bottom-Up):

- Vertical incision 1.5 to 2 cm long in the vaginal mucosa, starting 1.5 cm from the urethral meatus.
- Dissect laterally in the space between the vaginal wall and the urethral wall
- Make 2 suprapubic skin incisions approximately 1 cm long, 2 cm from the midline and 1 cm above the pubic symphysis.
- The loop of the pulling thread is inserted into the eye of the needle (Figure 2).

Figure 2: Instrument eye (a) and loop of the traction thread mounted onto the eye (b)


- Make the first passage of the needle vaginally, with a path along the posterior surface of the pubic symphysis.
- Take out the needle and collect the loop of the thread through the supra-pubic skin opening.
- Remove the needle respecting its path of insertion and its curvature.
- Repeat the same procedure for the other side.
- Place the tube close to the vaginal incision and remove the tube caps.
- Perform a follow-up cystoscopy.
- Make light traction movements to place the sling under the urethra.
- Tension adjustment: a good adjustment must allow the passage of the tip of a blunt scissor between sling and urethra.
- Perform a stress test.
- Section of the arms of the sling flush with the skin incisions.
- Closure of incisions with absorbable suture.

Trans-obturator approach (In-Out):

- Vertical incision 1.5 to 2 cm long in the vaginal mucosa, starting 1.5 cm from the urethral meatus.
- Dissection -with scissors up to the ischiopubic branch on each side.
- Make 2 skin incisions 1 cm long in the lower internal part of the obturator foramen at a horizontal level passing at least 1 cm over the urethral meatus and at most at the level of the clitoris.
- Introduction of the guide until contact with the ischiopubic branch.
- Passage from inside to outside of the needle, pulling thread mounted on the eye (Figure 2)..
- Recovery of the traction thread.
- Remove the needle respecting its path of insertion and its curvature.
- The operating time is reproduced on the contro lateral side.
- Place the tube close to the vaginal incision and remove the tube caps.
- Passage of the sling through the obturator foramen.
- Tension adjustment: a good adjustment must allow the passage of the tip of a blunt scissor between sling and urethra.
- Perform a stress test.
- Section of the arms of the sling flush with the skin incisions.
- Closure of incisions with absorbable suture.

Trans-obturator approach (Out-In):

- Vertical incision 1.5 to 2 cm long in the vaginal mucosa, starting 1.5 cm from the urethral meatus.
- Dissection with scissors to the ischiopubic branch on each side. The dissection should allow the index finger to be introduced to palpate the posterior part of the ischiopubic branch and the internal obturator muscle.
- Make 2 skin incisions 1 cm long in the lower internal part of the obturator foramen at a horizontal level passing at least 1 cm over the urethral meatus and at most at the level of the clitoris.
- Passage from outside to inside of the needle. The needle tip is localized with the index finger placed in the vaginal dissection.
- The loop of the pulling thread is inserted into the eye of the needle (Figure 2).
- The sling is exteriorized by pulling on the needle respecting its path of insertion and its curvature.
- The operating time is reproduced on the contro lateral side.
- Place the tube close to the vaginal incision and remove the tube caps.
- Tension adjustment: a good adjustment must allow the passage of the tip of a blunt scissor between sling and urethra.
- Perform a stress test.
- Section of the arms of the sling flush with the skin incisions.

- Closure of incisions with absorbable suture.

8 Side effects

Adverse events associated with mid-urethral sling placement can be de novo, persistent, worsened, transient or permanent. They can include but are not limited to:

- Partial or total failure of the intervention, leading to persistent incontinence or even worsening
- Erosion, extrusion or exposure of the sling through mucous membranes or surrounding organs (vaginal, urethral or bladder)
- Migration of the sling
- Infection of the sling or amplification of an existing urinary tract infection
- Urinary function dysfunction: overactive bladder, urgency, urinary obstruction / retention, pollakiuria, dysuria
- Sexual dysfunction: dyspareunia, spasm, loss of desire
- Local or systemic inflammation related to the presence of the sling
- Local irritation which may induce oedema or erythema
- Allergic reaction to polypropylene
- Vaginal discharge
- Discomfort, gene related to the presence of the device
- Transient or persistent pain

Surgery can also cause the following side effects:

- Perforation or laceration of surrounding vessels, nerves or organ
- Fistulas, dehiscence of the vaginal incision
- Bleeding, hematoma, hemorrhage, bruising, hematuria
- Scar tissue

In the presence of any of these side effects, it may be necessary to perform one or more corrective surgeries or even explant the prosthesis. These reoperations can induce undesirable effects of their own

9 Storage, handling and disposal

Store in the original packaging away from light and moisture.

Do not use the device beyond the expiration date shown on the label.

Return any defective product to the manufacturer or its distributor.

After use, dispose of the product and packaging in accordance with the regulations of the establishment, administration and / or local government.

Any explanted device should be returned to the manufacturer for analysis. In this case, do not clean or disinfect the device before shipment within a sealed waterproof packaging.

10 Information to be transmitted to the patient

The expected lifetime of the device is at least 15 years based on data accumulated on the device.

The patient should be handed the implant card and be informed by the surgeon of contraindications, possible side effects and post-surgical instructions to be followed as well as surgery failure risk if these ones are not followed.

SWING-BAND® sling is compatible with medical imaging techniques (radio, scanner and MRI etc.) it does not emit unintentional radiation, it is non-magnetic and radiolucent. The stainless steel of needles and guide are visible and compatible with X-rays but magnetic and therefore not compatible with MRIs

The patient must avoid any intense physical activity during the integration period of the prosthesis.

In the event of any undesirable side effect, the patient should contact her physician.

11 Implant card instruction

The implant card is available in the packaging. Before being handed to the patient, the implant card must be completed by the healthcare professional with the following information: patient name, surgery date and address of the healthcare institution and / or practitioner in page 2. See symbol signification in below Table.

THT BIO-SCIENCE

2 Cantignous
34220 Verreries-de-Moussans FRANCE
Tél. +33 (0)5 63 98 11 16 / Fax +33 (0)5 63 98 10 95
www.thtbio-science.com

SB IFU
Rev 01 25NOV2022
P 14 / 108


Figure 3 : Interior double-sided of the implant card



Figure 4 : Exterior double-sided of the implant card

Symbols	Signification
	CE mark and notified body identification number
	Manufacturer
	Country and date of manufacture
	Do not reuse
	Do not re-sterilize
	Do not use if package is damaged
	Expiry date
	Sterilized by ethylene oxide
	Catalog number
	Lot number
	Serial number
	Quantity

	Unique device identifier
	Device identifier
	Medical device
	Keep away from water
	Keep away from sunlight
	Unique sterile barrier system with inner protective packaging
	Consult the electronic instructions for use
	Patient identification
	healthcare institution or physician
	Date
2008	Date of affixing of the CE marking

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